
Thermoval®

duo scan

Gebrauchsanleitung
Mode d'emploi
Istruzioni per l'uso
Gebruiksaanwijzing
Instructions for use

- | | |
|---|--|
| <p>1 Mess-Sensor
Sonde de température
Sensore di misurazione
Meetsensor
Measuring sensor</p> <p>2 Aufsatz für Stirnmessung
Embout de mesure frontale
Cappuccio per misurazione sulla fronte
Kapje voor voorhoofdsmeting
Adapter for forehead measurement</p> <p>3 Thermometerspitze für Ohrmessung
Pointe de mesure auriculaire
Punta del termometro per misurazione nell'orecchio
Punt van de thermometer voor meting in het oor
Thermometer tip for ear measurement</p> <p>4 LED-Lampe zur Anzeige des Messvorgangs
Lampe LED pour visualiser l'affichage de la mesure
LED per indicazione del processo di misurazione
LED-lampje voor aanduiding van de meting
LED lamp to display the measuring process</p> | <p>5 SCAN-Knopf zum Start der Messung
Touche SCAN pour lancer la mesure
Pulsante SCAN per avviare la misurazione
Scan-knop voor het starten van de meting
SCAN button to start measurement</p> <p>6 Display
Ecran d'affichage
Display
Leesvenster
Display</p> <p>7 EIN-/AUS-Knopf zum Ein- und Ausschalten
Touche marche/arrêt pour allumer et éteindre l'appareil
Pulsante ON/OFF per accensione e spegnimento dell'apparecchio
Aan/uit-knop voor het aan- en uitschakelen
ON/OFF button for switching on or off the device</p> <p>8 Batteriefachdeckel
Couvercle du compartiment à piles
Coperchio del vano portabatterie
Deksel van het batterijkvakje
Cover of the battery compartment</p> |
|---|--|

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Paul Hartmannlaan 1
1480 Sint-Renelde/Saintes

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1. Preliminary remarks

 Please read these instructions for use carefully before first use as correct temperature measurement depends only on the appropriate use of the device. On account of external influences, wrong values may be measured otherwise. These instructions describe the individual steps of temperature measurement using Thermoal duo scan and contain important and helpful hints for the reliable determination of the body temperature. Be sure to keep these instructions for use for future reference.

Thermoal duo scan is a clinical thermometer, which uses infrared technology to measure the body temperature. Compared to conventional clinical thermometers, Thermoal duo scan has the advantage of producing the measurement faster. Moreover, the body temperature must no longer be measured orally, rectally or axillary.

By introducing the unit into the ear canal for measuring the temperature of the eardrum or by skin contact on the forehead, Thermoal duo scan can be used to reliably measure the body temperature. As the measuring process is shorter and the eardrum as well as the forehead is on the surface of the body, some measuring conditions have to be carefully observed in order to obtain a precise measurement.

The high measuring accuracy of Thermoal duo scan has been substantiated in clinical studies.

Table of contents

	Page
1. Preliminary remarks	47
2. Introduction	48
3. Benefits of the thermometer	48
4. General information on body temperature	48
5. Safety precautions	49
6. Initial operation of the device	49
7. Measurement of the body temperature	50
7.1. Ear measurement	51
7.2. Forehead measurement	52
8. Cleaning and maintenance of the device	54
9. Power supply and disposal notes	54
9.1. Batteries and disposal	54
9.2. Battery replacement	54
10. Error messages	54
11. Calibration check	55
12. Warranty conditions	55
13. Signs and symbols	55
14. Technical data	56
15. Legal requirements and guidelines	56
16. Accessories and spare parts	57
17. Contact information for customer queries	57

2. Introduction

Dear Customer, we are pleased that you have decided to purchase a clinical thermometer from HARTMANN. Thermoal duo scan is a quality product for measuring the human body temperature in the ear and on the forehead. It is ideally suited for measurements on children as from 6 months of age, but the thermometer can be used also on adults. Given correct application, the device ensures a fast and precise measurement of the body temperature in a very comfortable manner. We wish you all the best for your health.

3. Benefits of the thermometer

Two measuring methods: Ear and forehead measurement

Thermoal duo scan measures the temperature in the ear and on the forehead. Depending on your individual preference, you can choose between these two modes. The device recognizes automatically whether an ear measurement or forehead measurement is carried out.

High measuring accuracy ensured by the latest measuring sensor technology

The latest generation of the measuring sensors guarantees a very high measuring accuracy in ear and forehead measurements.

Quick measurement ensured by infrared technology

The measuring sensor picks up the infrared radiation emitted by the eardrum and/or the surface of the forehead and displays the measuring result within a second (three seconds in case of forehead measurement).

Easy to use

The device is very easy and safe to use having been provided with two buttons only. The O/I button is used to switch the device on (or off), and the SCAN button starts the measurement. In addition, the illustrations on the storage box explain how to carry out the measurements correctly.

Robust and reliable use ensured by high quality manufacturing

This device has been developed on the basis of the strict HARTMANN quality guidelines. Geared to the customer requirements, it fulfils the highest standards in stability.

4. General information on body temperature

The human body adjusts the body temperature to a setpoint value – it varies up to 1 °C throughout the day. In addition, the temperature in the body (core temperature) and the surface temperature of the skin show different temperatures. For this reason, there is no “normal” body temperature – it always depends on the measurement site. Here as well, the body temperature is influenced by the ambient temperature, age, stress, duration of sleep, hormones and physical activity.

Whilst glass thermometers and digital thermometers measure the temperature of the human body directly, thermometers using ear and forehead measurement determine the core temperature through the infrared radiation of the body. This may slightly deviate from the temperature obtained by rectal, oral or axillary measurement using digital thermometers, even when the measurement is carried out correctly. In the course of life, the average body temperature may decrease by up to 0.5° C.

Measurement table (in °C):

Designation	Measurement site		
	Ear / forehead	Anus	Mouth / armpit
Subnormal temperature	< 35.7	< 36.2	< 35.8
Normal temperature	35.8 – 36.9	36.3 – 37.4	35.9 – 37.0
Elevated temperature	37.0 – 37.5	37.5 – 38.0	37.1 – 37.5
Slight fever	37.6 – 38.0	38.1 – 38.5	37.6 – 38.0
Moderate fever	38.1 – 38.5	38.6 – 39.0	38.1 – 38.5
High fever	38.6 – 39.4	39.1 – 39.9	38.6 – 39.5
Very high fever	39.5 – 42.0	40.0 – 42.5	39.6 – 42.0

- The measurement in the ear and on the forehead can produce different results as the forehead temperature depends more on external influences than the temperature in the eardrum.
- The forehead measurement depends on external influence more than rectal, oral or axillary measurements.
- In case of several measurements after each other, slightly varying results are measured within the range of the measuring error tolerance as a rule.

5. Safety precautions

- This device may be used only for the measurement of human body temperature in the ear and on the forehead. The manufacturer shall not be held liable for any damage resulting from improper handling.
- Do not leave the device unattended near toddlers or persons who cannot operate it themselves.
- Do not expose the device to hard knocks or vibrations. Do not use the device, if you notice any kind of damage.
- Prior to every measurement, please check whether the glass lens of the sensor is clean, free from grease, and is not damaged.
- Please observe the storage and operating conditions mentioned in Chapter 14 "Technical data". Protect the device against dirt and dust, extreme temperatures, solar radiation and fluids.
- For cleaning, please follow the instructions mentioned in Chapter 8 "Cleaning and maintenance of the device".
- The device must not be altered, dismantled or repaired by the user.
- The use does not replace medical treatment in any way.
- If you inform your doctor about the temperature measured, please also mention that this temperature has been determined in the ear or on the forehead.

6. Initial operation of the device

Batteries are included in the scope of delivery and have been placed in the device already. However,

a contact strip interrupts the contact of the batteries to the device and thus guarantees the full battery capacity prior to the first use of the thermometer. Please carefully pull out the contact strip from the closed lid of the battery compartment. Subsequently the ThermoVal duo scan is ready for operation.

7. Measurement of the body temperature

ThermoVal duo scan measures the temperature on the forehead as well as in the ear. Both types of measurement can be used for the measurement of your own body temperature (self-measurement) as well as for the measurement of the body temperature of another person (foreign measurement).

1. Before starting the measurement, the device has to be switched on. For this reason, press the ON/OFF button (O/I). All symbols will be shown on the display for two seconds.
2. The device displays the temperature measured last for three seconds. At the same time you will see a symbol in the display indicating whether an ear measurement  or a forehead measurement  has been carried out.
3. Subsequently the temperature display will disappear and the thermometer symbol  will be faded in. Now the device is ready for measurement:

Important information to prevent inaccurate measured values:

- Please observe that the clinical thermometer and the person to be measured are in a room at


room temperature for at least 30 minutes prior to the measurement.

- In order to obtain the highest of accuracy in successive measurements, please wait for 30 seconds between two measurements.
- After sleeping, it is recommended to wait for a few minutes before measuring the temperature.
- Do not measure temperature directly after showering, swimming, et cetera, when the ears or the forehead are wet.
- Avoid eating, drinking or exercising prior to the measurement.
- Do not measure the temperature during or directly after breast feeding a baby.
- Remove the measuring device from the measuring area only after the end signal has been sounded and the LED lamp is no longer lit.
- In the following situations, we recommend to carry out three temperature measurements and to consider the highest value as measuring result:
 - Children below the age of three with a weak immune system (in particular, if the presence or absence of fever is significant).
 - In case of users, who are not familiar with the device, until they receive identical measured values.
 - In case of suspiciously low measured values.
- If the temperature measured is doubted and does not meet the feelings of the patient, it is recommended to repeat the measurement after some minutes. Please also refer to Chapter 8 "Cleaning and maintenance of the device". Furthermore, it is recommended to use another independent method of temperature measurement and/or to consult a doctor.

7.1 Ear measurement

The measurement can be carried out in the right-hand or left-hand ear. The thermometer measures the infrared rays emitted by the eardrum and the ear canal. These rays are picked up by the sensor and converted into temperature values. The best values are obtained when the infrared rays are picked up from the middle ear (eardrum). Rays which are emitted from surrounding tissue produce lower measured values.


If the cap for forehead measurement has been placed, please remove it by simply pulling it. For better grip, both sides of the cap have been provided with rough surfaces.


Stretch the ear canal by pulling the middle of the outer ear to the back and up at the same time. In case of children below the age of one year, pull the outer ear to the back only.

Carefully insert the thermometer tip into the ear canal until you feel a slight resistance.


Make sure that the tip is passed into the ear deep enough so that the ear symbol  appears.

Press the SCAN button to start measurement. At the beginning of measurement a brief tone sounds.

During the duration of the measurement, the LED over the SCAN button is lit. The ear measurement takes one second. When the measurement is finished, a long tone is sounded and the measuring result is displayed.

Please clean the tip of the thermometer and the glass lens of the sensor with the cleaning cloths enclosed and/or cotton buds.

The thermometer symbol  will appear three seconds after the measurement, and the device is ready for the next measurement.

Important information for ear measurement:


- Please make sure that the sensor is directed straight at the eardrum and not at the ear canal.
- If earwax accumulates on the measuring sensor, this may lead to inaccuracies of the temperature measurement or may even cause an infection between various users. For this reason, it is imperative that the measuring sensor is cleaned for every measurement. For notes on cleaning, please read the instructions in Chapter 8 "Cleaning and maintenance of the device".
- After cleaning the measuring sensor with alcohol-impregnated swabs, please wait for 15 minutes until carrying out the next measurement to permit the thermometer to reach the necessary operating temperature.
- The temperature values may vary from ear to ear. For this reason, please always use the same ear for measuring the temperature.
- No disposable protective caps are required for the measuring sensor. Please do not use disposable protective caps as they may distort the measuring result.
- If you have rested on an ear, the temperature in this ear may be elevated shortly after.
- The measurement must not be carried out in an ear which is suffering from inflammatory diseases as well as after ear injuries (such as eardrum injury) or during the healing phase after a surgical

procedure. If medication has been applied to an ear, this ear must not be used for measurement either. In these cases, please consult your attending doctor.


7.2 Forehead measurement

The thermometer measures the infrared radiation emitted by the skin in the forehead and temple area and the surrounding tissue. The radiation is picked up by the sensor and converted into temperature values. The most precise values are obtained when you scan the entire area starting from the middle of the forehead to the temple.

If the cap for forehead measurement has not been placed, please place it by simply pressing it on again.

Position the surface of the forehead adapter about 1 cm over the middle of the eyebrows. The thermometer has to have direct skin contact during the measuring process.

Press the SCAN button to start measurement. At the beginning of measurement a brief tone sounds.

Move the thermometer along the forehead towards the temple. During the duration of the measurement, the LED lamp over the SCAN button is lit. At the end of measurement, a long tone is sounded and the measuring result is displayed. The forehead measurement takes three seconds.

The thermometer symbol  will appear three seconds after the measurement, and the device is ready for the next measurement.

Important information for forehead measurement:

- We recommend measuring the temperature at the same position every time as otherwise the displayed values may vary.
- Remove hair, perspiration, cosmetics or dirt from the forehead to ensure an optimum measurement accuracy.
- Please observe that the forehead temperature may be influenced by staying outside for a longer period of time (in the cold in winter or in strong solar radiation in summer, for example) as well as by wearing a headgear.
- It is recommended that the forehead measurement is carried out after at least 30 minutes in room temperature.
- In order to ensure that the Arteria temporalis (temporal artery) is included during the measuring process, it is necessary to carry out the scan movement right up to the temple area, including the temple.

Switching over from ear to forehead measurement (or vice versa) is effected automatically by means of a contact sensor in the thermometer tip. The thermometer recognizes automatically which mode is operational.

8. Cleaning and maintenance of the device

The measuring sensor is the most important and most sensitive part of the device. In order to ensure most precise measurements, it always has to be clean and intact.

Clean the measuring tip and the sensor exclusively with the Soft-Zellin alcohol swabs supplied along with the device or with a soft moist cloth or with 70 percent alcohol (isopropanol). Residues from earwax can be removed carefully using cotton buds. Please do not use any aggressive cleaning agents or solvents as they may damage the device and could make the display intransparent. Please make sure that no fluid can penetrate the interior of the thermometer. To protect it from external influences keep the thermometer in the storage box.

9. Power supply and disposal notes

9.1 Batteries and disposal

- The two high-quality batteries included with the device guarantee at least 1,000 measurements. Use only high-quality batteries (see specification in Chapter 14 "Technical data").
- If you are using less efficient batteries, we can no longer guarantee 1,000 measurements.
- Never mix old and new batteries or batteries made by different manufacturers.
- Immediately remove exhausted batteries.
- If you do not intend to use the device for a longer period of time, you should remove the batteries to prevent possible leakage.
- Please help protect the environment: Do not dispose of batteries in household waste! Use designated collection points or municipal collection / recycling centres to dispose of hazardous

household waste.

This symbol on products and/or accompanying documents means that spent electronic products must not be mixed with regular household waste.

9.2 Battery replacement

This device has been provided with two alkaline batteries of LR03 (AAA) type. Please replace spent batteries by two new LR03 batteries as soon as the following battery symbol is indicated in the display:  .

Open the battery cover in the direction of the arrow on the rear side of the device by applying gentle pressure on the notch. Now insert the two "AAA" size batteries so that the positive (+) and negative (-) polarities match the "+" and "-" signs shown on the marking in the battery compartment. If the polarities do not match, the device will not work and the batteries may leak! Replace the lid of the battery compartment, applying gentle pressure until the hook snaps in place. The last measured value recorded in the memory will not be lost when changing the batteries.

10. Error messages

Thermoval duo scan is a clinically validated, high-quality product. Nonetheless error messages may occur, such as in the case of excessively high or low ambient temperature, if the measuring temperature is outside the range for human body temperatures, if the batteries are empty or in very rare cases, if a system error occurs. The description of errors is provided at the end of these instructions for use along with the explanation of the symbols on the display.

11. Calibration check

We recommend a calibration check at intervals of two years in the case of professionally used devices, for example in pharmacies, medical practices or hospitals. You should also observe the national regulations determined by the legislator. Calibration checks can be carried out either by the competent authorities or authorized maintenance providers against compensation.

12. Warranty conditions

We give a **3-year warranty** on this high-quality clinical ear and forehead thermometer from the date of purchase and in accordance with the following conditions.

Claims must be made during the warranty period. The date of purchase may be documented by the appropriately completed and stamped warranty certificate or proof of purchase.

Within the warranty period, HARTMANN shall replace or repair any faulty device components free of charge which were caused by material or manufacturing errors. This does not extend the warranty period.

This warranty is not applicable to damage caused by improper use or unauthorized interference. Parts that are subject to wear and tear (batteries, forehead adapter, storage box, packaging et cetera) are excluded from the warranty. Claims for compensation are limited to the value of the goods; compensation for subsequent damages is expressly excluded.

In warranty cases please send the device with cap, storage box and with the fully completed and

stamped warranty certificate direct, or via your dealer to the Customer Services department for your country.

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138, Australia

ZA – HARTMANN South Africa
2194 Johannesburg

13. Signs and symbols

Disposal note for electronic devices

Disposal note

Observe the operating instructions

Please note

Protection from electric shock (type BF)

14. Technical data

Type:	Infrared thermometer Thermoval duo scan
Measurement range:	Ear mode: 32.0 °C – 42.2 °C Forehead mode: 34.0 °C – 42.2 °C
Measuring accuracy:	Laboratory: ± 0.2 °C at 35.5 °C - 42.0 °C ± 0.3 °C at 32.0 / 34.0 °C – 35.4 °C and 42.1 °C – 42.2 °C in accordance with European Standard EN 12470-5
Display:	Liquid crystal display, 4 digits plus special symbols
Smallest display unit:	0.1 °C
Acoustic signal:	1 short signal tone when measuring process starts 1 long signal tone when measuring process ends
Optical signal:	Blue LED lamp is lit during the measuring process
Memory:	Last measured value
Background illumination:	The display is lit in blue when the device has been switched on and after measuring processes

Operating conditions:	Ambient temperature: 16 °C – 35 °C Relative humidity: 30 % – 85 %
Storage / transport conditions:	Ambient temperature: -25 °C – 55 °C Relative humidity: 30 % – 85 % Humidity for storage: 15 % – 85 %
Automatic switch-off:	About 1 minute after the last measurement
Power supply:	2 x 1.5 V Alkaline-manganese Mignon (AAA/LR03)
Battery capacity:	At least 1,000 measurements
Size:	150 mm (l) x 38 mm (w) x 40 mm (h)
Weight:	About 100 g (with battery)
Serial number:	See battery compartment

15. Legal requirements and guidelines

Thermoval duo scan complies with the requirements of the EC directive 93/42/EEC on medical devices (Medical Device Directive MDD) and bears the CE mark.

European Standard DIN EN 12470-5: Performance of infrared ear thermometers
ISO 80601-2-56

Portable and mobile high-frequency and communication devices can impair the functional capability of electronic medical devices.

16. Accessories and spare parts

Soft-Zellin alcohol swabs can be purchased in your pharmacy.
Code No. 999 979.

17. Contact information for customer queries

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138
Australia

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PAUL HARTMANN AG · 89522 Heidenheim, Germany
AU – PAUL HARTMANN Pty. Ltd. · Rhodes NSW 2138
BE – N.V. PAUL HARTMANN S.A. · 1480 Saintes/Sint-Renelde
CH – IVF HARTMANN AG · 8212 Neuhausen
GB – PAUL HARTMANN Ltd. · Heywood/Lancashire OL10 2TT
IT – PAUL HARTMANN S.p.A. · 37139 Verona