LANAFORM

# Skin Mass

SQUEEZING-ROLLING

SQUEEZING-ROLLING

TUNABLE INTENSITIES

PALPER ROULER

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PALPIER-ROLL-MASSAGE

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We would like to thank you for choosing "Skin Mass" by LANAFORM . This massaging device will prove to be an essential ally in regaining or keeping the figure of your dreams. "Skin Mass" not only helps eliminate orange peel skin but also enables you to regain skin that is smoother and more toned.

The positive effect of a Skin Mass massage lies in its technique, that is squeezing-rolling.

This massage technique, well known for many years to the institutes, guarantees better results when carried out over a long period. But as everyone knows, this practice can soon become long, tedious and costly.

With LANAFORM 's Skin Mass you can enjoy a massage when you want to depending on your schedule. Skin Mass has been specially developed for areas where grease collects rapidly such as the thighs, buttocks and stomach.

To start with, we recommend using Skin Mass and a slimming cream (not supplied) every day for about 10 minutes. When you first begin to see results, use it approximately three times per week.

## USING THE MASSAGER:

PLEASE READ ALL THE INSTRUCTION PRIOR TO USING THE "SKIN MASS", IN PARTICULAR THESE FEW BASIC SAFETY INSTRUCTIONS:

1. Only use this appliance for the purposes described in this manual.
2. This appliance is not intended for use by persons (including children) whose physical, sensorial or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they are supervised or instructed in the use of the appliance first by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
3. Please use the adaptor supplied with the appliance only. Using any other adapter risks harming the proper operation of the appliance.
4. Please do not use accessories which are not recommended by LANAFORM or which are not supplied with this appliance.
5. If it does not operate correctly, if it falls on the floor or if it is damaged or falls into water, have the appliance examined and repaired by the supplier or its after sales department.
6. Do not carry this appliance by its power supply lead or use this lead as a handle.
7. If the power lead is damaged, it must be replaced by a special or similar lead available from the supplier or after-sales department.
8. Do not use this product before going to bed. Massage has a stimulating effect which may inhibit sleep.
9. Never use this appliance in a room where aerosol products (sprays) are used or in a room where oxygen is administered.
10. Do not use the appliance on or close to the eyes or other sensitive areas.
11. Do not wrap the power lead around the appliance.
12. Make sure the product is not damaged before using it.
13. "Skin Mass" is a massaging appliance designed for non-professional use. Do not use this appliance as a substitute for proper medical care.
14. Keep the appliance in motion when using it.
15. The appliance is not watertight and cannot be used under the shower.

16. Use "Skin Mass" on clean, dry skin.
17. Never use the appliance if it has fallen into water.
18. Regular cleaning of the filter is essential because it can harm the suction.
19. Never operate the appliance under a blanket or pillow. Excessive heat may be created and cause a fire, electrocution or injury.
20. Always unplug the appliance from the power supply immediately after use and before cleaning.
21. Never insert objects into openings in the appliance.
22. Before unplugging the appliance, please be sure to set the control button to the "0" position.
23. An appliance should never be left unsupervised when it is plugged in.
24. If you are worried about your health, consult a doctor before using this appliance.
25. If you feel any pain while using this appliance, stop using it immediately and consult your doctor.
26. People fitted with a heart pacemaker should consult a doctor before using this device.
27. Never use this product directly on inflammation or swollen areas or where there is a skin rash.
28. Long use may lead to the production of excessive heat and shorten the life of the appliance. Consequently we advise using the appliance for a period of 10 minutes and then it should be left to cool down between sessions.

## DESCRIPTION OF THE "SKIN MASS":

### **A few words concerning the squeezing-rolling technique...**

Excess fats make skin look uneven and make the figure look displeasing, but the practice of squeezing and rolling helps disperse this fat and even prevents it from forming. It is a matter of carrying out a continuous action of pressure and pinching accompanied by progressive movement from bottom to top (in the direction of blood circulation) on the area to be treated, so as to cause a sort of rolling. The pressure and pinching knead and break the fat nodes and the heat created by this friction then helps the body remove them from the part being massaged. This technique also acts on blood microcirculation which promotes lymph drainage.

## HOW TO HOLD THE APPLIANCE

The "Skin Mass" operates via a 100~240 volts AC adapter. Connect the adapter to a plug socket and the other end to the female input connector located on the front of the appliance.

Hold the unit in the palm of your hand using the strap in the direction of the STOP button and the cable.

### ADJUSTING THE STRAP:

The strap can be adjusted to make using the massager more comfortable. Press the two ribbed buttons inside the strap and slide them away from the appliance. Raise the metal stem located inside the strap and select one of the three positioning slots. The strap can be adjusted on both sides of the appliance. After adjusting the strap, replace the clip sliding it into the appliance until you hear a click.

### OPERATION:

In order to understand how the product operates better, place it in front of you as shown in the photo.

1. "High" - "Low" wheel indicating the intensity.
2. ON / OFF (0) switch which also includes the 1, 2 and 3 functions of the squeezing-rolling system.
3. Indicator light to show which program is being used.

#### **Position 0 :**

The red light comes on when the appliance is plugged in and operational.

#### **Position 1 :**

This first position on your appliance starts the movement of the gentle roller that makes your skin more supple and gently raises it (The red light and green light 1 come on).

#### **Position 2 :**

This starts the adjustment of the "Skin Mass" suction mode as well as the movement of the gentle roller (The red light as well as the green lights 1 and 2 come on).

#### **Position 3 :**

This program adds vibration mode for a deep and even more pleasant massage (All the lights are on).

During operation of programs 2 and 3, you can adjust the suction intensity using the "high-low" wheel.

**NOTE: WE ALWAYS ADVISE STARTING YOUR SESSION USING THE LOWEST INTENSITY AND THEN PROGRESSIVELY INCREASING IT.**

**Stop Button:**

By pressing the central button you can switch off the suction and the movement of the rollers when you so wish.

**INSTRUCTIONS FOR USE:****Step 1**

Make sure that your skin is clean and dry. Apply cream to the areas to be treated, area by area so as to ensure that it does not penetrate into your skin before you can use the "Skin Mass".

Find a comfortable position and place the "Skin Mass" on your skin on the desired part of your body.

Switch the "Skin Mass" on and start your massage in position 1. The rollers will massage your skin and prepare it for the session.

Use the appliance from bottom to top towards your heart.

**Step 2**

Switch on the appliance's suction mode and adjust it. Always start with a low suction setting!

Now you can slowly adjust the intensity of the suction until you find the setting that is most comfortable for your massage.

The suction stop button must be used when you feel that suction has become too strong.

**Step 3**

After a few minutes use, switch on and adjust vibration mode for an even deeper and more effective massage.

**NOTE : WHEN USING THE SKIN MASS, ONLY THE SMALL ROLLER LOCATED UNDER THE DEVICE ROTATES. AUTOMATICALLY. ON THE OTHER HAND, THE LARGER ROLLER ROTATES WHEN THE SKIN MASS IS IN MOTION AND IN CONTACT WITH THE SKIN.**

**CLEANING AND STORAGE:**

- Use a dry cloth to wipe off any residues of cream.
- Always use a damp cloth to clean the appliance; do not use alcohol, acetone, abrasive products etc.

## ADVICE ON WASTE DISPOSAL

The packaging consists entirely of materials which do not pose any threat to the environment and which can be taken to your nearest sorting centre to be used as recycled material. Cardboard can be placed in a paper collection skip.

Packaging films should be delivered to your district sorting and recycling centre.

If you are no longer using the appliance, dispose of it in an environmentally friendly way and in accordance with the legal guidelines.

## LIMITED WARRANTY

LANAFORM guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM's control.

LANAFORM may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM's Service Centre. Any attempt to maintain this product by a person other than LANAFORM's Service Centre will render this guarantee void.

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