

SANITAS

Pèse-
personne



SBF 14

SBF 14

Pèse-personne en verre

Plateforme carrée en verre de sécurité

Boutons à effleurement modernes

Ecran LCD facile à lire: taille des chiffres 26 mm

Electrodes en acier brossé inoxydable

Affiche le poids, le taux de graisse, l'eau corporelle, la masse musculaire et la masse osseuse

Graduation:

- 0,1% pour la graisse et l'eau corporelles et le taux de muscle
- 100 g pour le poids

5 niveaux d'activité

10 mémoires utilisateurs avec reconnaissance automatique

Sélecteur kg, lb, st

Capacité de charge: 180 kg

Indication en cas de surcharge

Arrêt automatique

Démarrage: fonction "Step-On" (en montant dessus)

1 pile x 3V CR2032

Dimensions: 30 x 30 x 2,3 cm

Poids: ca. 1550 g

Garantie 3 ans

UE: 4 / Transportkarton: -

EAN-Nr.: 4211125 73591 3

Art.Nr.: 735.91

SBF 14

Glass diagnostic scale

Square platform of safety glass

Modern sensor buttons

Easy to read LCD-display: digit size 26 mm

Electrodes of brushed stainless steel

Display of body weight, body fat, body water, muscle ratio and bone mass

Graduation:

- 0.1% for body fat, body water and muscle ratio
- 100 g for body weight and bone mass

5 activity levels

12 user memories with automatic user recognition

Change over to kg, lb, st

Capacity: 180 kg

Overload indication

Automatic switch-off

Switch-on technology: Step-On

1 x 3V CR2032 battery included

Product measurements: 30 x 30 x 2.3 cm

Product weight: approx. 1550 g

3 year guarantee

Sales unit: 4 / Shipping unit: -

EAN no.: 4211125 73591 3

Item no.: 735.91