



german|engineering

BC 40 Tensiomètre poignet

Mesure entièrement automatique de la tension artérielle et du pouls au niveau du poignet

Moyenne de toutes les mesures enregistrées et de la tension du matin et du soir des 7 derniers jours

1 x 60 Emplacements de mémoire

Indicateur OMS* (échelle des valeurs selon les directives de l'OMS*)

Détection de l'arythmie: met en garde en cas de turbulence du rythme cardiaque

Produit médical

Manchette pour une taille de poignet de: 12,5 - 21,5 cm

Date et heure / Arrêt automatique

Message en cas d'erreur d'utilisation

Mesure automatique par inflation: pour une mesure rapide et confortable

Passeport de tension artérielle

Indicateur du niveau des piles

Pochette de rangement pratique

2 piles AAA de 1,5V incluses

Dimensions (L/H): 76 x 70 x 23 mm

Garantie 3 ans

UE: 4

Colisage: 4 x 4

N°EAN: 42 11125 65061 2

N°d'art.: 650.61

BC 40 Blood pressure monitor

Fully automatic blood pressure and pulse measurement on wrist

Average of all saved measured values and morning and evening blood pressure for the last 7 days

1 x 60 memory spaces

With WHO* indicator (classification of values according to WHO* guidelines)

Arrhythmia detection: warning in case of possible heart rhythm disturbance

Medical product

Cuff size for wrist circumference 12.5 - 21.5 cm

Date and time / automatic switch-off

Signal in case of application error

Automatic measurement during inflation: for fast and comfortable measurement

Blood pressure passport

Low battery indicator

Practical storage box

Inkl. 2 x 1.5V AAA batteries

Product measurement (L/W/H): 76 x 70 x 23 mm

3 years guarantee

Sales Unit: 4

Shipping carton: 4 x 4

EAN no.: 42 11125 65061 2

Item no.: 650.61

*WHO / *OMS = World Health Organization (Organisation Mondiale de la Santé)